




# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT



# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT





In this graph you can see that there are two color shaded areas of data. These areas are referred to as plots of information from your continuous glucose monitoring (CGM) device. The blue plot is your pump and sensor information from the dates in Date Range A. Because this is the most recent information downloaded from your pump, an average sensor glucose (SG) line is calculated and shown as a dotted black line in the middle. The dark shading in blue represents 25-75% of all your sensor readings, meaning this is where most of your glucose readings have been. Remember, your CGM records up to 288 SG values on a daily basis, from those 288 values, 25-75% of them are represented in the darker shade. The remaining or excess data are in the 0-90% range shown within the solid blue line, outside of the darker shaded area.


Your data from Date Range B, is colored in orange behind the blue plot. This section of the report should be reviewed with your HCP to see progress from your last visit or your last device settings change. Do you see less shading in the blue plot below 70mg/dL compared to the orange plot? This is a good discussion to start with your HCP to see if you are having high frequency with low glucose.

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT




## Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)



<b>3</b>	<b>A</b>	13.0	12.0	12.0
	<b>B</b>	13.0	12.0	12.0

**A** Hypoglycemic  
**1** 8:15 PM (1 occurrence)



X

For easier comparison and to review with your healthcare provider, your carb ratio settings from both Date Range A and B have been displayed.

Hypoglycemic patterns (1) # Episodes (per day): 0.1  
 10:50 PM - 1:00 AM

	<b>A</b>	<b>B</b>
<b>Statistics</b>		
SmartGuard (per week)	100% (7d 00h)	100% (7d 00h)
Manual Mode (per week)	0% (00h)	0% (00h)
Sensor Wear (per week)	96% (6d 18h)	96% (6d 17h)
Average SG ± SD	155 ± 41 mg/dL	156 ± 46 mg/dL
GMI***	7.0%	7.0%
Coefficient of Variation (%)	26.5%	29.8%
Low / High SG Alerts (per day)	0.0 / 0.1	0.5 / 6.7
Average BG	159 ± 52 mg/dL	156 ± 54 mg/dL
BG / Calibration (per day)	2.6 / 2.5	3.6 / 3.1
Total daily dose (per day)	34 units	40 units
Bolus amount (per day)	17U (48%)	13U (33%)
Auto Correction amount (per day)	8U (47%)	13U (97%)
Auto Basal / Basal amount (per day)	18U (52%)	27U (67%)
Set Change	Every 3.0 days	Every 2.7 days
Reservoir Change	Every 3.5 days	Every 4.7 days
Meal (per day)	7.5	0.4
Carbs entered (per day)	107 ± 17 g	104 ± 31 g
Active Insulin time	2:00 hrs	2:00 hrs

\*\*\* Glucose Management Indicator


This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT



## Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)




Carb Ratio (g/U)

**A** 13.0 **B** 13.0

**4** Hypoglycemic patterns (1) # Episodes (per day): 0.1

**1** 8:15 PM - 8:50 PM (1 occurrences)



### SmartGuard Exits **A**

No Calibration	0
SmartGuard max delivery	0
SmartGuard min delivery	0
BG required for SmartGuard	0
Sensor Algorithm Underread	0
Sensor Updating	0
No SG values	0
Sensor Expired	0
SmartGuard disabled by user	0
Prolonged Suspend	0
SmartGuard Warm Up	0
Unidentified	0

Note hypoglycemic (low glucose) patterns numbered during the times they occurred. Each occurrence is an episode that lasted 30 minutes in duration. These patterns apply to Date Range A only.

Average BG	159 ± 52 mg/dL	156 ± 54 mg/dL
BG / Calibration (per day)	2.6 / 2.5	3.6 / 3.1
Total daily dose (per day)	34 units	40 units
Bolus amount (per day)	17U (48%)	13U (33%)
Auto Correction amount (per day)	8U (47%)	13U (97%)
Auto Basal / Basal amount (per day)	18U (52%)	27U (67%)
Set Change	Every 3.0 days	Every 2.7 days
Reservoir Change	Every 3.5 days	Every 4.7 days
Meal (per day)	7.5	0.4
Carbs entered (per day)	107 ± 17 g	104 ± 31 g
Active Insulin time	2:00 hrs	2:00 hrs

\*\*\* Glucose Management Indicator


This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT




## Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)



Carb Ratio (g/U)	13.0	12.0	12.0
<b>A</b>	13.0	12.0	12.0
<b>B</b>	13.0	12.0	12.0

Hypoglycemic patterns (1)	# Episodes (per day):	0.1	Hyperglycemic patterns (1)	# Episodes (per day):	0.1
<b>A</b>			<b>5</b>		
<b>1</b> 8:15 PM - 8:50 PM (1 occurrences)			<b>2</b> 10:50 PM - 1:00 AM		



SmartGuard Exits	<b>A</b>	<b>B</b>
No Calibration	0	0
SmartGuard max delivery	0	0
SmartGuard min delivery	0	0
BG required for SmartGuard	0	0
Sensor Algorithm Underread	0	0
Sensor Updating	0	0
No SG values	0	0
Sensor Expired	0	2
SmartGuard disabled by user	0	0
Prolonged Suspend	0	0
SmartGuard Warm Up	0	0
Unidentified	0	0

Note hyperglycemic (high glucose) patterns numbered during the times they occurred. Each occurrence is an episode that lasted 30 minutes in duration. These patterns apply to Date Range A only.

Set Change	Every 3.0 days	Every 2.7 days
Reservoir Change	Every 3.5 days	Every 4.7 days
Meal (per day)	7.5	0.4
Carbs entered (per day)	107 ± 17 g	104 ± 31 g
Active Insulin time	2:00 hrs	2:00 hrs

\*\*\* Glucose Management Indicator


This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT




## Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)



Carb Ratio (g/U)	12 AM - 6 PM	6 PM - 12 AM	12 AM - 6 PM
<b>A</b>	13.0	12.0	12.0
<b>B</b>	13.0	12.0	12.0

Hypoglycemic patterns (1)	# Episodes (per day):	Hyperglycemic patterns (1)	# Episodes (per day):
<b>A</b>	0.1	0.1	0.1
<b>1</b> 8:15 PM - 8:50 PM (1 occurrences)		<b>2</b> 10:50 PM - 1:00 AM	



\*\*\* Glucose Management Indicator


This part of the report shows time in range. Time in Range is the amount of time spent in target; for this representation, the green range (70-180mg/dL). You have the ability to compare Date Range A and Date Range B as an evaluation of your glucose management with your HCP. Any percentage in the less optimal ranges of 40-69mg/dL or 181-400mg/dL should be discussed with your HCP for fine tuning and adjustment.

<b>A</b>	<b>B</b>
100% (7d 00h)	100% (7d 00h)
0% (00h)	0% (00h)
96% (6d 18h)	96% (6d 17h)
155 ± 41 mg/dL	156 ± 46 mg/dL
7.0%	7.0%
26.5%	29.8%
0.0 / 0.1	0.5 / 6.7
159 ± 52 mg/dL	156 ± 54 mg/dL
2.6 / 2.5	3.6 / 3.1
34 units	40 units
17U (48%)	13U (33%)
8U (47%)	13U (97%)
18U (52%)	27U (67%)
Every 3.0 days	Every 2.7 days
Every 3.5 days	Every 4.7 days
7.5	0.4
107 ± 17 g	104 ± 31 g
2:00 hrs	2:00 hrs

This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center


# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT

Medtronic Assessment and Progress  
 (A) 12-02-2019 - 12-15-2019 (14 Days) (B) 12-16-2019 - 12-29-2019 (14 Days)



Carb Ratio (g/U)  
 (A) 13.0  
 (B) 13.0

Hypoglycemic patterns (1) # Episodes (per week)  
 (A) 1 8:15 PM - 8:50 PM (1 occurrences)  
 (B) 1 8:15 PM - 8:50 PM (1 occurrences)



\*\*\* Glucose Management Indicator

This table can be used to help you understand the frequency and causes of your SmartGuard™ exits if you are wearing the MiniMed 780G system. Review this section with your HCP to uncover behaviors that will assist you in avoiding some of these SmartGuard™ exits.

SmartGuard Exits	(A)	(B)
No Calibration	0	0
SmartGuard max delivery	0	0
SmartGuard min delivery	0	0
BG required for SmartGuard	0	0
Sensor Algorithm Underread	0	0
Sensor Updating	0	0
No SG values	0	0
Sensor Expired	0	2
SmartGuard disabled by user	0	0
Prolonged Suspend	0	0
SmartGuard Warm Up	0	0
Unidentified	0	0


Statistics	(A)	(B)
SmartGuard (per week)	100% (7d 00h)	100% (7d 00h)
Manual Mode (per week)	0% (00h)	0% (00h)
Sensor Wear (per week)	96% (6d 18h)	96% (6d 17h)
Average SG ± SD	155 ± 41 mg/dL	156 ± 46 mg/dL
GMI***	7.0%	7.0%
Coefficient of Variation (%)	26.5%	29.8%
Low / High SG Alerts (per day)	0.0 / 0.1	0.5 / 6.7
Average BG	159 ± 52 mg/dL	156 ± 54 mg/dL
BG / Calibration (per day)	2.6 / 2.5	3.6 / 3.1
Total daily dose (per day)	34 units	40 units
Bolus amount (per day)	17U (48%)	13U (33%)
Auto Correction amount (per day)	8U (47%)	13U (97%)
Auto Basal / Basal amount (per day)	18U (52%)	27U (67%)
Set Change	Every 3.0 days	Every 2.7 days
Reservoir Change	Every 3.5 days	Every 4.7 days
Meal (per day)	7.5	0.4
Carbs entered (per day)	107 ± 17 g	104 ± 31 g
Active Insulin time	2:00 hrs	2:00 hrs

This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center



# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT

Medtronic




- Use these statistics to quickly glance at your progress between both date ranges. Percentage of SmartGuard™ and Manual Mode are shown in days and hours and should be reviewed with your HCP to optimize your glucose management.
- Take a look at your sensor wear and your average SG. This is an average of all your sensor glucose readings. Does your average SG appear to be in the range of your glucose control goals?

8 Statistics		<b>A</b>	<b>B</b>
SmartGuard (per week)		100% (7d 00h)	100% (7d 00h)
Manual Mode (per week)		0% (00h)	0% (00h)
Sensor Wear (per week)		96% (6d 18h)	96% (6d 17h)
Average SG ± SD		155 ± 41 mg/dL	156 ± 46 mg/dL
GMI***		7.0%	7.0%
Coefficient of Variation (%)		26.5%	29.8%
Low / High SG Alerts (per day)		0.0 / 0.1	0.5 / 6.7
Average BG		159 ± 52 mg/dL	156 ± 54 mg/dL
BG / Calibration (per day)		2.6 / 2.5	3.6 / 3.1
Total daily dose (per day)		34 units	40 units
Bolus amount (per day)		17U (48%)	13U (33%)
Auto Correction amount (per day)		8U (47%)	13U (97%)
Auto Basal / Basal amount (per day)		18U (52%)	27U (67%)
Set Change		Every 3.0 days	Every 2.7 days
Reservoir Change		Every 3.5 days	Every 4.7 days
Meal (per day)		7.5	0.4
Carbs entered (per day)		107 ± 17 g	104 ± 31 g
Active Insulin time		2:00 hrs	2:00 hrs

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT

Medtronic




- The Glucose Management Indicator (GMI) is calculated based on mean sensor glucose and reflects glycemic control for that reporting period. GMI provides a useful measure for connecting CGM metrics to laboratory A1C and reinforces your engagement with your HCP to a personalized diabetes management plan. Reports with at least 14 days in the reporting period and 10 days of sensor wear, are needed to generate this value.
- Coefficient of Variation (%CV) is a measurement of glycemic variability similar to standard deviation (SD). Talk with your healthcare provider about relative risk of lows if %CV is 36% or greater.<sup>1</sup>
- Your Low/High SG alerts (per day) explain how many sensor alerts you are experiencing on average during the reporting period. This can also help assist you and your HCP if some alarms might need adjustment in settings or can be turned off.

<sup>1</sup> Diabetes Care Aug 2019, 42 (8) 1593-1603; DOI: 10.2337/dci19-0028

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT

Medtronic



- Did your average blood glucose (BG) in Date Range A improve from Date Range B?
- Are you maintaining 3 to 4 BG readings/day according to the American Diabetes Association (ADA) guidelines?
- Do you calibrate 3-4 times/day for optimal sensor performance or enter a blood glucose (BG) as the system has requested to remain in SmartGuard™?

# MINIMED™ 780G SYSTEM | ASSESSMENT & PROGRESS REPORT




## Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)

Percentile comparison

25-75% 0-90%

Average **A**




12.0  
12.0  
# Episodes (per day): 0.1

**A** **B**

Week)	100% (7d 00h)	100% (7d 00h)
Week)	0% (00h)	0% (00h)
Week)	96% (6d 18h)	96% (6d 17h)
	155 ± 41 mg/dL	156 ± 46 mg/dL
	7.0%	7.0%
tion (%)	26.5%	29.8%
ts (per day)	0.0 / 0.1	0.5 / 6.7
	159 ± 52 mg/dL	156 ± 54 mg/dL
er day)	2.6 / 2.5	3.6 / 3.1

- This table shows you the distribution of your insulin usage between bolus and basal insulin. Take a look at the insulin total daily dose. How much insulin do you use on average per day? You can use this number to see how much insulin is needed on a monthly basis. Auto Correction amount and percentage is calculated from Bolus amount (per day).
- Set Change and Reservoir Change can help you focus on the frequency of your infusion set and reservoir changes. Are you following the recommended routine given by your HCP and to the labeling of your specified infusion set and reservoir?



\*\*\* Glucose Management Indicator

Sensor Expired	0	2 ●●
SmartGuard disabled by user	0	0
Prolonged Suspend	0	0
SmartGuard Warm Up	0	0
Unidentified	0	0

<b>11</b>	Total daily dose (per day)	34 units	40 units
	Bolus amount (per day)	17U (48%)	13U (33%)
	Auto Correction amount (per day)	8U (47%)	13U (97%)
	Auto Basal / Basal amount (per day)	18U (52%)	27U (67%)
	Set Change	Every 3.0 days	Every 2.7 days
	Reservoir Change	Every 3.5 days	Every 4.7 days
	Meal (per day)	7.5	0.4
	Carbs entered (per day)	107 ± 17 g	104 ± 31 g
	Active Insulin time	2:00 hrs	2:00 hrs


This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center



### Assessment and Progress

**A** 12-02-2019 - 12-15-2019 (14 Days) **B** 10-11-2019 - 10-24-2019 (14 Days)

Percentile comparison



- Use this section to monitor how many carbs per day you are eating. Everyone is different so check with your HCP or contact a dietitian for appropriate nutritional recommendations for daily carbohydrate consumption. Remember to enter all your carbs into the pump. Snacks are also included in your carbs entered calculation. Your meals (per day) will also include those snacks.
- Active insulin time is the amount of time it takes for food or correction insulin to lower your blood glucose. On average, active insulin time is set to 3 hours, however, confirm with your HCP if this amount of time is appropriate for you as everyone metabolizes insulin differently.

# Episodes (per day): 0.1

	<b>A</b>	<b>B</b>
Week)	100% (7d 00h)	100% (7d 00h)
Week)	0% (00h)	0% (00h)
Week)	96% (6d 18h)	96% (6d 17h)
	155 ± 41 mg/dL	156 ± 46 mg/dL
	7.0%	7.0%
tion (%)	26.5%	29.8%
ts (per day)	0.0 / 0.1	0.5 / 6.7
	159 ± 52 mg/dL	156 ± 54 mg/dL
er day)	2.6 / 2.5	3.6 / 3.1
er day)	34 units	40 units
day)	17U (48%)	13U (33%)
ount (per day)	8U (47%)	13U (97%)
amount (per day)	18U (52%)	27U (67%)
	Every 3.0 days	Every 2.7 days
	Every 3.5 days	Every 4.7 days

<b>12</b> Meal (per day)	7.5	0.4
Carbs entered (per day)	107 ± 17 g	104 ± 31 g
Active Insulin time	2:00 hrs	2:00 hrs

\*\*\* Glucose Management Indicator

This report is compatible with the Ambulatory Glucose Profile calculations used by the International Diabetes Center

## CARELINK™ SOFTWARE IMPORTANT SAFETY INFORMATION

The CareLink™ software is intended for use as a tool to help manage diabetes. The purpose of the software is to take information transmitted from insulin pumps, glucose meters and continuous glucose monitoring systems, and turn it into CareLink™ reports. The reports provide the information that can be used to identify trends and track daily activities – such as carbohydrates consumed, meal times, insulin delivery, and glucose readings. NOTE: CareLink™ report data is intended for use as an adjunct in the management of diabetes only and NOT intended to be relied upon by itself. Patients should consult their healthcare providers familiar with the management of diabetes prior to making changes in treatment. For more details, please consult <https://www.medtronicdiabetes.com/important-safety-information> and the appropriate CareLink™ User Guide at <https://www.medtronicdiabetes.com/download-library>.

# Medtronic

18000 Devonshire Street  
Northridge, CA 91325  
USA  
800-646-4633

Toll-free: 800-328-2518  
24-Hour Technical Support for  
physicians and medical professionals)

950M23911-011 25NOV2020 Medtronic. All rights reserved.  
Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.  
All other brands are trademarks of a Medtronic company.

[www.medtronicdiabetes.com](https://www.medtronicdiabetes.com)