Medtronic

ENLITE™

Glucose Sensor • Capteur de glucose
Glukose-Sensor • Sensor de glucosa
Glucosesensor • Sensore del glucosio
Glukosesensor • Glukoosisensori
Glukossensor • Glukosesensor
Sensor de glicose • Senzor za glukozo
Senzor pro sledování glykemie • Glükózszenzor
Sensor glukozy • Glikoz Sensörü
Glukózoný senzor • Αισθητήρας γλυκόζης
Сенсор глюкозы
グルコースセンサー • 葡萄糖传感器
חישן סוכר • جهاز استشعار الغلوكوز

Indications for use

The Enlite glucose sensor (sensor) is intended for use with Medtronic Diabetes (Medtronic) glucose sensing systems to continuously monitor glucose levels in persons with diabetes. Refer to the Enlite Serter User Guide for instructions on inserting the sensor. The sensor has a maximum life of 144 hours (six days). The 144-hour life span of the sensor begins when the insulin pump or monitor receives the first METER BG NOW alert.

Contraindications

None known.

Warnings

The Enlite Serter (Serter) does not work the same as other Medtronic insertion devices. **The Serter injects the sensor into the insertion site when the button is released**, not when the button is pressed. Be sure to read the entire Enlite Serter User Guide before attempting to insert the sensor. Failure to follow directions may result in pain or injury.

The sensor may create special needs regarding your medical conditions or medications. Bleeding, swelling, irritation or infection at the insertion site are possible risks associated with inserting the sensor and sometimes result from improper insertion and maintenance of insertion site. Discuss these conditions and medications with your doctor before using the sensor.

Remove the sensor if you experience unexplained fevers, or develop inflammation, redness, soreness or tenderness at the insertion site. Check the insertion site often for these conditions and to ensure that the sensor is still in place. Consult your healthcare professional if these conditions persist.

This device contains an automatically retracting needle which may cause minimal blood splatter when the needle is removed from the sensor. When removing the needle from someone other than yourself, do so at arms length and pointed away from self and others. Keep device within sight at all times to avoid accidental needle stick injury.

Sensors are sterile and nonpyrogenic, unless the package has been opened or damaged. Always inspect the packaging for damage prior to use. Do not use any of the sensors if the sterile package has been opened or damaged.

Watch for bleeding at the insertion site. If bleeding occurs, do not attach the transmitter or recorder to the sensor. Apply steady pressure using a sterile gauze or clean cloth for up to 3 minutes.

If bleeding stops, attach the transmitter or recorder to the sensor.

If bleeding does *NOT* stop, DO NOT attach the transmitter or recorder to the sensor.

1. Remove the sensor and discard.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation and treat accordingly.
3. Insert a new sensor in a different location.

This product contains small parts and may pose a choking hazard for young children.

Precautions

Do not reuse sensors. Reuse of the sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation or infection.

Always wash your hands with soap and water before opening the sensor package and handling the sensor.

Healthcare professionals and caregivers should use universal precautions when handling the sensor. Avoid touching any sensor surfaces that will come in contact with the body—that is, the sensor, needle, adhesive surfaces and tape.

Do not insert the sensor through tape and only use alcohol to clean the insertion site. Do not use any other type of skin preps prior to insertion.

Establish a rotation schedule for choosing each new sensor site. Dispose of used sensors and introducer needles in a sharps container after a single use. Do not clean or resterilize, and do not try to extract the needle from the needle housing.


Where to insert the sensor

Refer to the Enlite Serter User Guide for instructions on inserting the sensor.

WARNING: Failure to follow directions may result in pain or injury.

CAUTION: Make sure the sensor site is at least 7.5 cm (3 inches) away from the insulin pump insertion site or manual injection site. When replacing the sensor, select a new site that is at least 5 cm (2 inches) from the previous site.

Choose an insertion site that has an adequate amount of subcutaneous fat. Shown here are the best body areas (shaded) for sensor insertion.

Do not insert the sensor in muscle or areas constrained by clothing or accessories; areas with tough skin or scar tissue; sites subjected to rigorous movement during exercise; in sites on the belt or waistline; or within 5 cm (2 inches) of the navel.

Reagents

The sensor contains two biological reagents: glucose oxidase and human serum albumin (HSA). Glucose oxidase is derived from *Aspergillus Niger*, purified and dried according to Type VII-S guidelines and cross-linked to HSA, purified and dried albumin fraction V, derived from pasteurized human serum with gluteraldehyde. Since less than 0.4 µg of glucose oxidase and less than 0.7 µg of HSA is used to manufacture each sensor, the risks of tissue reactions and viral transmission are considered minimal.

Storage and handling

CAUTION: Do not freeze the sensor. Sensors should not be stored in direct sunlight, extreme temperatures, or humidity.

Store sensors at room temperature between +36°F to +86°F (+2°C to +30°C). If the sensor is stored in a cool environment, allow the sensor to warm to room temperature to prevent condensation. Do not store sensors at temperatures lower than +36°F (+2°C).

Discard sensor after “Use by” on label, or if the package is damaged or the seal is broken.

	Use by
	Do not reuse
	Follow instructions for use.
REF	Catalogue number
LOT	Batch code
STERILE R	Sterilized using irradiation
	Temperature limitation
	Date of manufacture
	Open here
	Manufacturer

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©2010, Medtronic Minimed, Inc. all rights reserved. Enlite™ は Medtronic MiniMed, Inc. の商標です。
Medtronic Ltd.
Medtronic Österreich GmbH
Medtronic B.V.