# Guardian<sup>™</sup> Connect

Application User Guide

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Australia: Medtronic Australasia Pty. Ltd. Tel: 1800 668 670

Azerbaijan: Isomed Tel: +994 (12) 464 11 30

Bangladesh: Sonargaon Healthcare Pvt Ltd. Mobile: (+91)-9903995417 or (+880)-1714217131

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**Malaysia:** Medtronic International Ltd. Tel: +603 7946 9000

Middle East and North Africa: Regional Office Tel: +961-1-370 670

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Thailand: Medtronic (Thailand) Ltd. Tel: +662 232 7400

Türkiye: Medtronic Medikal Teknoloji Ticaret Ltd. Sirketi. Tel: +90 216 4694330

**USA:** Medtronic Diabetes Global Headquarters 24 Hour HelpLine: +1-800-646-4633 To order supplies: +1-800-843-6687

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## Introduction

Thank you for choosing Medtronic as your diabetes management partner.

The Guardian<sup>™</sup> Connect application (app) is a component of the Guardian Connect Continuous Glucose Monitoring (CGM) system. The Guardian Connect CGM system helps you manage your diabetes by:

- Recording your glucose values throughout the day and night
- Showing the effects that diet, exercise, and medication can have on your glucose levels
- Giving you additional tools, such as alerts, the ability to record diet and exercise, and access to CareLink<sup>™</sup> Personal

This user guide is designed to help you understand the setup and operation of your Guardian Connect app.

## System description

CGM is a sensor glucose monitoring tool that uses a glucose sensor placed below your skin to continuously measure the amount of glucose in your interstitial fluid. It allows users to track patterns in glucose concentrations and to identify episodes of low and high glucose. The Guardian Connect CGM system includes the **Guardian Connect app (CSS7200)**, an app that can be downloaded onto a mobile device with a compatible operating system version. A list of compatible operating system versions is available in the app store found on your mobile device. The app displays sensor glucose data, and also provides a user interface for sensor calibration, for entering data such as exercise and meals, and for uploading information to the CareLink Personal website.

The Guardian Connect app displays alerts if a glucose level reaches, falls below, or rises above set values. The app also provides the interface for manually entering fingerstick blood glucose calibration values obtained from a commercial blood glucose meter. The app then stores the sensor glucose data so that it can be analyzed to track patterns.

The Guardian Connect CGM system also includes the **Guardian Connect transmitter (MMT-7821)** and the **Enlite<sup>™</sup> sensor (MMT-7008)**. This guide describes how to set up and use the Guardian Connect app. For details on using the other components of the Guardian Connect CGM system, please see the corresponding device user guides.

## Intended use

The Guardian Connect app is intended for continuous or periodic monitoring of glucose levels in the interstitial fluid under the skin, in persons with diabetes mellitus.

The Guardian Connect app is intended for use with a compatible consumer mobile electronic device. It allows users to track patterns in glucose concentrations and to possibly identify episodes of low and high glucose. The Guardian Connect app displays alerts if a glucose level reaches, falls below, or rises above set values.

Sensor glucose values displayed on the screen are not intended to be used directly for making therapy adjustments, but rather to provide an indication of when a meter blood glucose measurement may be required.

## Contraindications

None known.

## Safety Warnings

Caution: Any changes or modifications to devices not expressly approved by Medtronic could interfere with your ability to operate the equipment and void your warranty.

The Guardian Connect app requires that you have a functioning mobile electronic device with the correct settings (such as a supported operating system and Bluetooth® enabled) for accurate operation. A non-functioning mobile device, or incorrect settings on your mobile device, may prevent the Guardian Connect app from properly displaying alerts.

## Precautions

You must test your blood glucose levels at least two times per day, or as indicated by the system. If the Guardian Connect app indicates that your sensor glucose is not within your glucose target range, check your blood glucose using your blood glucose meter.

Taking medications with acetaminophen while wearing the sensor may falsely raise your sensor glucose readings. The level of inaccuracy depends on the amount of acetaminophen active in your body and may be different for each person. Always use blood glucose meter readings to verify your glucose level before making therapy decisions.

## Assistance

Please contact your local help line or representative for assistance. Refer to the Medtronic Diabetes International Contacts list at the beginning of this user guide for the help line or representative in your area.

## How to use this guide

The following table describes terms and conventions used in this guide.

Conven- tion	Description		
Toggle	Indicates that the same feature on the screen can be used to switch between two op- tions. For example, "Toggle an alert on" means that you slide a switch right to turn on an alert. To turn it off, you need to slide the same switch left.		
Bold	Indicates an item on the screen that you select with your finger or tap to open.		
>	A shorthand to indicate a series of selections you make on the screen. For example, Alert Settings > Rate Alerts means that you need to tap Alert Settings, and then on the next screen tap Rate Alerts.		
Note	Provides additional helpful information.		
CAUTION	Notifies you of a potential hazard which, if not avoided, may result in minor or moderate injury or damage to the equipment.		
WARNING	Notifies you of a potential hazard which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.		

## **Guardian Connect app setup**

When you open Guardian Connect for the first time, the app walks you through the setup process. Simply follow the instructions on your screen. For more information, see your welcome card.

## Pairing your transmitter

Follow the instructions on your screen to pair your transmitter to your mobile device.

#### New sensor setup

Follow the instructions in the user guide for inserting the sensor. Then follow the onscreen instructions to complete the sensor start up.

Note: You must connect your transmitter to your sensor before completing set up.

For details on reconnecting an existing sensor, see *Reconnecting the existing sensor, on page 22.* 

#### Inserting the sensor

Refer to your serter user guide for instructions on how to insert the sensor.

#### Connecting the transmitter to the sensor

Refer to your transmitter user guide for instructions on how to connect the transmitter to the sensor.

## Home screen

The Home screen appears when you open the Guardian Connect app on your mobile device.

**Note:** When you open your Guardian Connect app for the first time, there will be no sensor information displayed on the Home screen. Your first sensor glucose reading appears after you have successfully paired your transmitter and calibrated your sensor.

Item Description		
Calibration	Displays the Calibration screen where you enter a blood glucose (BG) meter read- ing for sensor calibration. For details, see <i>Calibrating your sensor, on page 8</i> .	
	The color of the ring and the size of its displayed segment indicate calibration sta- tus and the approximate time left until your next sensor calibration is due. When your sensor is fully calibrated, the ring is solid green. As the time for your next sen- sor calibration approaches, segments of the ring disappear, and the color of the ring changes, until a red blood drop appears instead of the ring, as shown below:	
	O (0) (0) (0) (0) (0) (0) (0) (0) (0) (0)	
	A full green circle indicates 12 hours remaining to the next calibration. Yellow indi- cates three hours remaining. Red indicates one hour remaining. When the red blood drop appears, an alert is displayed instructing you to calibrate your system. The question mark indicates that the calibration status is unknown. Three dots indi- cate that a calibration currently is not permitted.	
Rate of change arrows	Display your glucose trend and the rate at which the most recent sensor glucose level has risen or fallen. For details on setting the rate of change for falling and rising sensor glucose values, see <i>Setting rate alerts, on page 14</i> .	
Menu screen	Provides access to the menu screen, which displays the Guardian Connect Sys- tem status information and allows you to set up your sensor, define alert settings, view the logbook, and access CareLink Personal settings. For details, see <i>Menu</i> <i>screen, on page 6</i> .	
Continuous sen- sor glucose trace	Displays your current and previous SG readings. Tap any point on the graph to view the details of the selected sensor glucose reading or event in a box that appears above. For more navigating tips, see <i>Navigating through the sensor glucose graph, on page 15.</i>	
System status icons	Display the corresponding status icon if the transmitter battery is empty, the life of your sensor has expired, the connection with the transmitter is lost, or the status of these features is unknown. For details, see <i>System status icons, on page 7</i> .	
Time interval	Displays preset time intervals of 3, 6, 12, and 24 hours. Tap this icon to switch be- tween time intervals.	
Current sensor glucose reading	Displays your current sensor glucose reading, which the transmitter calculates and sends wirelessly to the Guardian Connect app.	
High and low sensor glucose alert limits	Display a line showing the high and low sensor glucose limits on the sensor graph. The orange line indicates your sensor glucose high limit; the red line indicates your sensor glucose low limit.	
Sensor status message	Displays the latest active status notification. If a current sensor glucose reading is not available, the sensor status message will appear where the sensor glucose information would be displayed.	

ltem	Description	
Sensor glucose graph	Displays a blue line representing your sensor glucose readings over a selected time interval. It also indicates your high and low glucose limits. For details on the sensor graph, see <i>Sensor graph, on page 15</i> .	
Day indicator	Displays the date of the data shown in the graph. As you scroll further back into historical data, the date is displayed as Today, Yesterday, days of the previous week.	
Add event	Displays the events screen where you can enter additional information, such as exercise, blood glucose readings, meals you eat, or insulin you take. Certain blood glucose readings entered here may be also used for sensor calibration. For details, see <i>Events, on page 17</i> .	

#### Menu screen

The Menu screen appears when you tap 😑 on the Guardian Connect Home screen.

Item	Description	
Close (X)	Closes the menu screen and returns you to the Home screen.	

ltem	Description		
System status bar	Displays icons that show a quick status of your Guardian Connect system, includ- ing transmitter battery, sensor life, and transmitter communication. For details on the individual icons and their status, see <i>System status icons, on page 7</i> .		
Logbook	Displays a history of alerts and events that have occurred in the Guardian Connect System. For details, see <i>Logbook, on page 21</i> .		
Alert Settings	Display the options for setting glucose monitoring alerts and calibration reminders. For details, see <i>Alert settings, on page 11</i> .		
Sync to Care- Link	Provides access to a menu for configuring data uploads to CareLink Personal, for changing your user login information, and for granting your care partner(s) access to your diabetes information via the CareLink Personal website.		
Sensor Setup	Displays the option for pairing a transmitter to your mobile device. Also allows for unit selection (glucose and carbs) where available, per country regulations.		
About	Provides the user with software version information, access to the user guide, and reference to the end user license agreement.		

#### System status icons

The system status icons appear at the top of the Menu screen. These icons provide a way for you to quickly check the status of your system. If any condition becomes critical and requires your immediate attention, the corresponding status icon also is displayed on the Guardian Connect Home screen next to the Calibration icon.

Icon name	Description		
Transmitter battery	The charge level of your transmitter battery. The color and the fill level of the icon indicate the status. When your battery is full, the icon is solid green. The question mark indicates that the battery status is unknown. As the battery is used, the icon changes, as shown below:		
Sensor life	The number of days remaining in the life of your sensor. The number and the color of the icon indicate the status. When you insert a new sensor, the icon is solid green. As your sensor life is used, the number on the icon indicates the number of days left before you need to replace your sensor, and the color of the icon changes as shown below:		
	The question mark indicates that the sensor life is unknown. When there is less than one day left, the 1-day icon remains red, and the text below it indicates the number of hours remaining.		

Icon name	Description	
Transmitter communication	The communication status between your transmitter and the Guardian Connect app:	
	$\bigotimes$ — the transmitter is active and connected	
	A method with the transmitter is not paired to your mobile device, or Bluetooth is off	

## Calibrating your sensor

Calibration is the process of entering a blood glucose meter reading to calculate sensor glucose values. You must calibrate your sensor regularly to ensure you continue to receive accurate sensor glucose data. For details, see *Calibrating guidelines, on page 9*.

**Note:** The Guardian Connect app requires up to a two-hour warm-up, from the time the sensor is connected to the transmitter, before you can calibrate your sensor. The calibration icon remains hidden until the sensor has completed its warm-up.

Whenever a sensor calibration is due, you will receive an alert instructing you to calibrate your sensor now. After you calibrate the sensor, it may take up to five minutes to receive updated sensor glucose values.

Required calibration	When?	Description
First	Within two hours after connecting a new sensor.	Sensor will take up to two hours to become operational. As soon as it is initialized, you will be asked to calibrate your sensor.
Second	Within six hours after your first calibration.	If you do not calibrate for more than six hours after the first calibration, your transmitter stops calculating your sensor glucose values until you calibrate your sensor.
Subsequent	Within 12 hours after your second calibration and at least every 12 hours thereafter.	If you do not calibrate for more than 12 hours, your transmitter stops calculating your sensor glucose values until you calibrate your sensor. You must calibrate the sensor at least every 12 hours throughout its life. For details, see <i>Calibrating guidelines, on page 9</i> .

The following table describes when a sensor calibration is required.

**Note:** You may also receive additional calibration alert to let you know that another calibration is required to improve sensor performance. If you do not calibrate after receiving an additional calibration alert, your transmitter stops calculating your sensor glucose values until you calibrate your sensor.

## Calibrating guidelines

Follow these guidelines for best sensor calibration results:

- Calibrate three to four times at regular intervals throughout the day to improve accuracy. If necessary, set a calibration reminder as explained in *Setting calibration reminders, on page 15.*
- Enter your blood glucose meter reading into the app immediately after testing your blood glucose.
- Always wash your hands before testing your blood glucose.
- Use only your fingertips when obtaining blood samples for calibration.
- Avoid using an old blood glucose reading or reusing blood glucose readings from previous calibrations.
- If for some reason calibration is unsuccessful, you must wait at least 15 minutes before attempting another calibration.

**Note:** If your blood glucose meter readings are significantly different from your sensor glucose readings, wash your hands and calibrate again.

#### Entering blood glucose meter readings for calibration

You can enter a blood glucose meter reading for calibration from the Calibration screen as described in this section or when you enter a blood glucose meter reading on the Events screen, as described in *Entering blood glucose meter readings, on page 17*.

To enter a blood glucose meter reading for sensor calibration:

- 1 Take a blood glucose meter reading.
- 2 On the Guardian Connect Home screen, tap the calibration icon.
- 3 Using the number pad, enter the blood glucose value between 2.2 mmol/L and 22.2 mmol/L (40 mg/dL and 400 mg/dL).
- 5 Tap **Calibrate** at the top. Confirm the value, and then tap **Calibrate sensor XXX mmol/L (mg/dL)**.

The app returns to the Home screen, and a blood drop appears on the sensor glucose graph at the time of the entry. Your sensor glucose reading appears within five minutes after calibration.

## Alerts

An alert makes you aware of a situation that may require your immediate attention. You must always resolve the underlying problem that triggered an alert. On your Guardian Connect app, Guardian Connect alerts appear as colored popup messages. When the Guardian Connect app is not open, or your device is locked, the style and appearance of alerts depends on the notification settings on your mobile device.

- CAUTION: Do not turn off notifications or enable the Do Not Disturb setting on your mobile device. Doing so may result in missing important alerts requiring your immediate attention.
- CAUTION: Do not force the Guardian Connect app to quit. Doing so will stop transmission of sensor glucose data, and you will not receive any alerts. If the app is running in the background on your mobile device, you will still receive sensor glucose data and may still receive alerts depending on your mobile device's notification settings.
- CAUTION: Medtronic cannot guarantee that your mobile device will always allow for the Guardian Connect app to run in the background. If the mobile device automatically closes your Guardian Connect app, you may receive a 'Lost Communication' notification. Particularly when running several applications at once, you should not rely on alerts from Guardian Connect and you should periodically check if the app is still running.

Depending on the kind of underlying problems, Guardian Connect alerts may be grouped into glucose monitoring alerts (displayed as red and orange popup messages within the app) and error and system status alerts (displayed as blue popup messages).

**Error and system status alerts** notify you when a new sensor is required, communication with the transmitter is lost, calibration is required, or other conditions exist that require your attention to ensure successful operation of the Guardian Connect system. Except for calibration reminders (Calibrate By alert), you do not need to set these alerts, as the system automatically generates them when the related conditions occur.

**Glucose monitoring alerts** notify you if your sensor glucose values are changing at a particular rate, are approaching, or have reached or exceeded a low or high limit that you have specified. To receive these alerts, you need to set your glucose limits and define the alert settings as described in *Alert settings, on page 11*.

The following graph illustrates the different high and low glucose alerts you can use.

## Alert settings

The following table describes the types of glucose alerts and other settings available on the Alert Settings screen of the Guardian Connect app:

Setting	Description		
High Alerts	Use High Alerts to define when you will be notified if your sensor glucose is approaching or has reached the high limit. For details, see <i>Setting high alerts, on page 11</i> .		
Low Alerts	Use Low Alerts to define when you will be notified if your sensor glucose is approaching or has reached the low limit. For details, see <i>Setting low alerts, on page 13</i> .		
Rate Alerts	Use Rate Alerts to define when you will be notified if your sensor glucose is rising or falling at a specified rate. For details, see Setting rate alerts, on page 14.		
Snooze Time	Use Snooze Time to define when you will be reminded that an alert condition still ex- ists. After an alert is received and snoozed, you will be alerted again only if the alert condition still exists after the snooze time you have set. For details, see <i>Setting alert</i> <i>snooze time, on page 15</i> .		
Calibration Reminder	Use Calibration Reminder to define when you will be reminded that a sensor calibration is due. For details, see <i>Setting calibration reminders, on page 15</i> .		

#### Setting high alerts

The following table describes the high alert settings:

High alert setting	Description	
High Limit	Your high limit is the value on which your other high settings are based. Your high limit can be set from 5.6 mmol/L to 22.2 mmol/L (100 mg/dL to 400 mg/dL). On your sense glucose graph, your high limit appears as an orange horizontal line at the value that y set.	
Alert on High	When Alert on High is on, your system displays a High Sensor Glucose alert when your sensor glucose value reaches or exceeds your high limit.	
Alert Be- fore High	When Alert Before High is on, you will receive a High Predicted alert any time the sen- sor glucose is predicted to reach your high limit. This makes you aware of potential hig glucose levels before they occur.	
Time Be- fore High	The Time Before High option is available only when the Alert Before High feature is o This option determines when you will receive a High Predicted alert, if your sensor glu cose values continue to increase at the current rate of change. You can set a time be tween ten minutes and one hour.	

You can define the same high alert settings for the entire day or use different settings during the day and at night.

To set your high alerts:

- 1 On the Guardian Connect Home screen, tap and select **Alert Settings** > **High Alerts**.
- 2 Toggle **All Day** to switch between all-day alerts and different alerts for day and night. If you choose to set different alerts for day and night, two sets of settings appear.
- 3 For daytime alerts, tap **Day starts at** and set the desired start time for the selected day period.
- 4 Tap High Limit and set the desired limit between 5.6 mmol/L and 22.2 mmol/L (100 mg/dL to 400 mg/dL).
  This tagges Alert on Using an

This toggles Alert on High on.

- 5 Return to the High Alerts screen and tap **Save**.
- 6 If you want to receive alerts when your sensor glucose is approaching your high limit:
  - a. Toggle Alert Before High on.
  - b. Tap **Time Before High** and set the desired time when you want to receive a High Predicted alert.
- 7 If you are setting different alerts for day and night, tap **Night starts at**, select the time, and repeat steps 4 and 5 to set your nighttime alerts.
- 8 Tap **Save** at the top of the screen to save your settings. Then tap **Alert Settings** to define the other alert settings, such as snooze time for the alerts you have just set.

9 When you have completed setting up your alerts, tap **Home** to return to the Guardian Connect Home screen.

#### Setting low alerts

The following table describes the low alert settings:

Low alert setting	Description	
Low Limit	Your low limit is the value on which your other low settings are based. Your low limit can be set from 2.8 mmol/L to 5 mmol/L (50 mg/dL to 90 mg/dL). On your sensor glucose graph, your low limit appears as a red horizontal line at the value that you set.	
Alert on Low	When Alert on Low is on, your system displays a Low Sensor Glucose alert when your sensor glucose value reaches or falls below your low limit.	
Alert Be- fore Low	When Alert Before Low is on, you will receive a Low Predicted alert any time the sensor glucose is predicted to reach your low limit. This makes you aware of potential lows before they occur.	
Time Be- fore Low	Time Before Low is available only when Alert Before Low is on. Time Before Low deter- mines when you will receive a Low Predicted Alert, if your sensor glucose values contin- ue to decrease at the current rate of change. You can set a time between ten minutes and one hour.	

You can define your low alert settings for the entire day or use different settings during the day and at night.

To set your low alerts:

- 1 On the Guardian Connect Home screen, tap = and select **Alert Settings** > **Low Alerts**.
- 2 Toggle **All Day** to switch between all-day alerts and different alerts for day and night.

If you choose to set different alerts for day and night, two sets of settings appear.

- 3 For daytime alerts, tap **Day starts at** and set the start time for the selected day period.
- 4 Tap Low Limit and set the desired limit between 2.8 mmol/L and 5 mmol/L (50 mg/dL to 90 mg/dL). This toggles Alert on Low on.
- 5 Return to the Low Alerts screen and tap Save.
- 6 If you want to receive alerts when your sensor glucose is approaching your low limit:

- a. Toggle Alert Before Low on.
- b. Tap **Time Before Low** and set the desired time when you want to receive a Low Predicted alert.
- 7 If you are setting different alerts for day and night, tap **Night starts at**, select the time, and repeat steps 4 and 5 to set your nighttime alerts.
- 8 At the top of the screen, tap **Alert Settings** to save your settings and to define the other alert settings, for example, snooze time for the alerts you have just set.
- 9 When you have completed setting up your alerts, tap **Home** to return to the Guardian Connect Home screen.

#### Setting rate alerts

Rate alerts notify you when your sensor glucose is rising (Rise Alert) or falling (Fall Alert) equal to or faster than a specified rate. These alerts help you understand how your glucose levels are affected, for example, by meals or exercise.

On the Guardian Connect Home screen, these rapidly rising or falling trends are indicated by arrows, as shown in *Home screen, on page 4*. The more arrows, the higher the rate of change.

To set your rate alerts:

- 1 On the Guardian Connect Home screen, tap and select **Alert Settings** > **Rate Alerts**.
- 2 Toggle Rise Alert on.

Three options appear, with the number of arrows corresponding to a specific rise rate that is equal to or faster than the rate you have selected to trigger the Rise Alert. One arrow signals a .056 mmol/L/minute (1 mg/dL/minute) rise; two arrows signal a .111 mmol/L/minute (2 mg/dL/minute) rise; and three arrows signal a .167 mmol/L/minute (3 mg/dL/minute) rise.

- 3 Tap the arrow option with the rise rate you want to use.
- 4 Toggle Fall Alert on.

Three options appear, with the number of arrows corresponding to a specific fall rate that is equal to or faster than the rate you have selected to trigger the Fall Alert. One arrow signals a .056 mmol/L/minute (1 mg/dL/minute) fall; two arrows signal a .111 mmol/L/minute (2 mg/dL/minute) fall; and three arrows signal a .167 mmol/L/minute (3 mg/dL/minute) fall.

- 5 Tap the arrow option with the fall rate you want to use.
- 6 At the top of the screen, tap **Save** > **Alert Settings** to return to the Alert Settings screen. Then tap **Home** to return to the Guardian Connect Home screen.

#### Setting alert snooze time

If you have set your high, low, or rate alerts and want to be reminded if the condition that had triggered the alerts still exists after a certain period of time, you can set a snooze time for your alerts. You can set snooze time for high and rise alerts that is different from snooze time for low and fall alerts.

To set your alert snooze time:

- 1 On the Guardian Connect Home screen, tap and select **Alert Settings** > **Snooze Time**.
- 2 Tap **High and Rise Alerts** and set the desired amount of time before you are reminded of the existing alert condition.
- 3 Tap **Low and Fall Alerts** and set the desired amount of time before you are reminded of the existing alert condition.
- 4 At the top of the screen, tap **Save** > **Alert Settings** to return to the Alert Settings screen. Then tap **Home** to return to the Guardian Connect Home screen.

#### Setting calibration reminders

You can set a reminder for yourself that sensor calibration is due within a specified period of time. For example, if you set your reminder to 30 minutes, you will receive a Calibrate By alert 30 minutes before you need to enter a BG meter reading for calibration.

To set a calibration reminder:

- 1 On the Guardian Connect Home screen, tap and select **Alert Settings** > **Calibration Reminder**.
- 2 Toggle **Reminder** on.
- 3 Set the desired amount of time to be alerted before sensor calibration is due.
- 4 At the top of the screen, tap **Save** > **Alert Settings** to return to the Alert Settings screen. Then tap **Home** to return to the Guardian Connect Home screen.

## Sensor graph

The sensor graph displays your current sensor glucose reading, which the transmitter sends via Bluetooth to your Guardian Connect app. It also allows you to view a history of sensor glucose readings and events you have entered.

#### Navigating through the sensor glucose graph

Use the following tips to navigate the Guardian Connect sensor graph data:

- Swipe the center of the graph right and left to view historical data and to return to the current graph location.
- Pinch and stretch the center of the graph with your fingers to zoom in and out on the graph data.
- Double-tap the graph to view the graph at the predefined 3-hour, 6-hour, 12-hour, and 24-hour zoom levels.
- Tap the horizontal time axis to center the selected time point on the graph and to view the details of the selected sensor glucose reading or event in an information box that appears above the graph. For details, see *Graph information boxes, on page 16*.
- Tap the vertical sensor glucose values axis twice to return to the current sensor glucose reading displayed on the graph.

#### Graph information boxes

When you tap any point on the graph, the tapped time is marked by a vertical cursor on the graph, and a box with the information about the tapped event or sensor glucose reading appears above the graph. The following figure provides an example:

For sensor glucose readings, information boxes display the sensor glucose value, its date and time, and rate of change arrows, similarly to the current sensor glucose value information displayed on the Home screen. If there is no sensor glucose value for the selected time point, the information box displays the sensor status message for that point. For event markers, information boxes display the details specific to each event type. For details, see *Events, on page 17*.

You can drag the cursor across the graph, skipping at five-minute intervals to pinpoint a specific sensor glucose value on the graph.

Information boxes appear for a few seconds and then close.

## Events

Events help you capture information that may affect your glucose levels. You can use the Events screen on the Guardian Connect app to enter and save certain types of events.

Event icon	Event name	Description
$\bigcirc$	Blood Glucose	Blood glucose meter readings. These can be used either to calibrate the system or simply to manage your diabetes without calibrating the system.
6	Insulin	The type and amount of insulin you have delivered.
	Meal	The amount of carbohydrates you eat or drink.
S	Exercise	The intensity and duration of your exercise routine.
ß	Other	This event can be used to enter any other information relevant to your diabetes management. For example, you can record information such as when you take medications, feel ill, or are under stress.

## **Note:** Make it a practice to enter events when they happen. If you make a mistake when entering an event, you can always delete the event and enter it again.

#### Entering blood glucose meter readings

If you measure your blood glucose, for example, when you eat or when your blood glucose is rising or falling rapidly, you can enter these measurements into the Guardian Connect app.

You also have the option of using the entered blood glucose meter reading for calibration purposes, if calibration is allowed at the time when you entered the event. To use the blood glucose meter reading for calibration, the value must be between 2.2 mmol/L and 22.2 mmol/L (40 mg/dL and 400 mg/dL). When calibrating your sensor from the Blood Glucose screen, you may enter your current blood glucose value, or values that are up to 12 minutes old.

**Note:** You can enter a blood glucose value between 1.1 mmol/L and 33.3 mmol/L (20 mg/dL and 600 mg/dL) into the Events feature. For blood glucose calibration, you can enter a value between 2.2 mmol/L and 22.2 mmol/L (40 mg/dL and 400 mg/dL).

To enter a blood glucose meter reading into the Events feature:

- 1 Measure your blood glucose with a blood glucose meter.
- 2 Tap + at the bottom of the Guardian Connect Home screen.
- 3 Tap on the Events screen.
- 4 If you need to change the date or time for the entry, tap **Time** and make the appropriate changes.
- 5 Using the number pad, enter the blood glucose meter reading.
- 6 Make sure the value displayed above the number pad is correct. If the value is incorrect, tap 🕲 to clear it and enter the correct value.
- 7 Tap **Save** at the top of the screen.
- 8 If you receive a message instructing you to calibrate now:
  - To update the sensor with the entered blood glucose meter reading, select Calibrate sensor XXX mmol/L (mg/dL).
  - If you do not want to update the sensor with this reading, select **Just save**.
  - If you need to make a correction or need to cancel, tap **Cancel**.

The app returns to the Home screen, and **o** appears on the graph at the selected time.

#### Entering insulin injection information

If you have delivered insulin using an insulin pump, insulin pen, or a syringe, you can enter the amount of insulin you injected into the Guardian Connect app.

To enter the type and amount of insulin you have delivered:

- 1 Make a note of the amount and type of injection you have taken.
- 2 At the bottom of the Guardian Connect Home screen, tap +.

- 3 On the Events screen, tap (6).
- 4 If you need to change the date or time for the entry, tap **Time** and make the appropriate changes.
- 5 Using the number pad, enter the insulin amount.
- 6 Make sure the value displayed above the number pad is correct. If not, tap 🗷 to clear it and enter the correct value.
- 7 Tap **Type** and select the type of insulin you have taken.
  - Rapid acting
  - Long acting
- 8 Tap **Done** at the top.

The app returns to the Home screen, and appears on the graph at the selected time.

#### Entering meal information

You can use the Guardian Connect app to record information about the carbohydrates you eat or drink with meals or snacks.

To enter your meal information:

- 1 Determine the total amount (in grams) of carbohydrates in the meal, snack, or drink that you plan to consume.
- At the bottom of the Guardian Connect Home screen, tap +.
- 3 On the Events screen, tap ().
- 4 If you need to change the date or time for the entry, tap **Time** and make the appropriate changes.
- 5 Using the number pad, enter the amount of carbohydrates you have consumed.
- 6 Make sure the value displayed above the number pad is correct. If not, tap 🗷 to clear it and enter the correct value.
- 7 Tap **Done** at the top.

The app returns to the Home screen, and nappears on the graph at the selected time.

#### Entering exercise information

You can use the Guardian Connect app to enter information about your exercise regimen. Make sure you are consistent and enter the marker either before or after each time you exercise.

To enter your exercise information:

- 1 Make a note about how long you exercised (duration), and how difficult or easy the exercise was (intensity).
- 2 At the bottom of the Guardian Connect Home screen, tap +.
- On the Events screen, tap
- 4 If you need to change the date or time, tap **Time** and set the desired day and time for this exercise entry.
- 5 Tap **Duration** and set the time you have spent exercising. The entered value appears above the number pad.
- 6 Tap **Intensity** and select Low, Medium, or High to indicate how intensely you have exercised.
- 7 Tap **Done** at the top.

The app returns to the Home screen, and spears on the graph at the selected time.

#### **Entering other events**

You can use the Guardian Connect app to enter events other than blood glucose measurements, insulin injections, carbohydrates consumed, and exercise information. For example, you can record information such as when you take medications, feel ill, or are under stress.

To enter other events:

- 1 At the bottom of the Guardian Connect Home screen, tap +.
- 2 On the Events screen, tap 0.
- 3 Using the text field, enter the relevant information.
- 4 Tap **Done** at the top.

The app returns to the Home screen, and papears on the graph at the selected time.

## Logbook

The Logbook screen displays a history of alerts and events that occurred on the selected day, with the most recent entries at the top of the list.

To view Logbook entries:

- 1 On the Guardian Connect Home screen, tap **=** and select **Logbook**.
- 2 Do any of the following to view the desired information:
  - Tap **Alerts** or **Events** to filter the list by the specific type. You also can select **All** to view the entire list.
  - Swipe down and up on the list to view the entries.
  - If you wish to delete an event entry, swipe it left and tap Delete.

Note: You cannot delete alerts and calibration events.

- Tap the desired entry to expand it and view the details. You can also use and varrows at the top of the expanded view to scroll through the list entries.
- 3 If you are viewing a specific event or alert on a details screen, tap **Logbook** > **Home** to return to the Home screen. If you are on the main Logbook screen, tap **Home** to return to the Home screen.

## Syncing your data to CareLink Personal

Guardian Connect allows you to sync your data to CareLink Personal automatically or upload information for CareLink reports as needed.

This automatic Sync to CareLink feature sends data displayed in the Guardian Connect app to the CareLink Connect tab of the CareLink Personal website approximately every five minutes when an Internet connection to CareLink is available. This feature also automatically sends sensor history information used to create CareLink reports roughly every 24 hours. This information can be viewed by you or a care partner on the CareLink Personal website at carelink.minimed.eu. Please note that your mobile device must be connected to the Internet to send data to CareLink. If using a cellular connection, your provider's data rates may apply.

The Sync to CareLink feature defaults to on when you set up your Guardian Connect app. If the Sync to CareLink toggle is turned off, the app will no longer send sensor information to CareLink Personal.

The **Upload Now** button allows you to immediately send sensor history data to CareLink Personal for generating reports.

## Sending a care partner access to CareLink

Now that your Guardian Connect app is synced with CareLink, you can invite a family member, friend, or care partner to track your diabetes information on the CareLink website. Visit carelink.minimed.eu on your personal computer or tap the **Manage Care Partners** button within the app to access the CareLink Personal website. There you can create a username and temporary password for each care partner. Your care partners use this login information to access your diabetes information on the CareLink Personal website.

**Note:** Your care partners will use the CareLink login information you create for them to access the CareLink website.

For more information on the CareLink Connect tab, please refer to the *CareLink Connect User Guide* found in the CareLink Connect tab of CareLink Personal.

## Disconnecting the transmitter from the sensor

Always refer to your transmitter user guide for instructions on disconnecting the transmitter from the sensor.

## Reconnecting the existing sensor

If necessary, you can reconnect your transmitter to an existing sensor. Simply connect your transmitter to your sensor. When the app detects the connection, confirm that the sensor is an existing sensor. It may take a few seconds to establish a connection when connecting an existing sensor. If you reconnect an existing sensor, the sensor will go through another warm-up period before you can calibrate.

## Icon table

	Manufacturer
(((••)))	Radio communication
REF	Catalogue or model number
8	Follow instructions for use.
	Bluetooth <sup>®</sup> wireless technology or Bluetooth <sup>®</sup> enabled
(MR)	Magnetic Resonance (MR) Unsafe: keep away from magnets and strong electromagnetic fields.

	RF device distributed in Australia
EC REP	Authorized representative in the European community
<b>C € 0459</b>	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC (NB 0459).
CONF	Configuration or unique version identifier

©2016, Medtronic MiniMed, Inc. All rights reserved.

MiniMed<sup>™</sup>, Guardian<sup>™</sup>, Enlite<sup>™</sup>, and CareLink<sup>™</sup> are trademarks of Medtronic MiniMed, Inc.

Bluetooth® is a registered trademark of Bluetooth SIG, Inc.

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©2016, Medtronic MiniMed, Inc. Tutti i diritti riservati.
Bluetooth® è un marchio registrato di Bluetooth SIG, Inc.
Medtronic B.V. Earl Bakkenstraat 10 6422 PJ Heerlen The Netherlands