- (en) Application Sites
- (es) Zonas de aplicación
- (de) Applikationsstellen
- (it) Posizioni di applicazione
- (fr) Sites d'application
- (nl) Toepassingszones
- (da) Påsættelsessteder
- (fi) Käyttökohdat
- (no) Brukssteder
- (sv) Applikationsställen
- (pl) Miejsca stosowania
- (cs) Místa aplikace
- (el) Σημεία εφαρμογής
- (pt) Locais de Aplicação
- (sl) Mesta aplikacije
- (hu) Felhelyezési területek

Unomedical
A Convatec Company

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Indications for Use

The i-port Advance™ injection port is indicated for patients who administer, or to whom is administered, multiple daily subcutaneous injections of physician prescribed medications, including insulin. The device may remain in place for up to 72 hours to accommodate multiple drug injections without the discomfort of additional needle sticks. The i-port Advance™ may be used on a wide range of patients, including adults and children. Model numbers for the device are MMT-100 and MMT-101.

Product Description

The i-port Advance™ is a subcutaneous injection port, which combines the injection port with an aid for insertion. The injection port is an integrated part of the inserter. Two injection port models are available: a 6mm or 9mm catheter length. The i-port Advance™ is delivered assembled and ready for use as a single, sterile unit.

Contraindications

- DO NOT use the same i-port Advance™ for longer than 72 hours.
- DO NOT attempt to reapply the i-port Advance™ after use.
- The i-port Advance™ is neither intended nor indicated for use with blood or blood products.
- DO NOT use i-port Advance™ for any purpose other than the indications specified by your healthcare provider.
- DO NOT use a needle longer than 5/16 in. (8mm) and /or thicker than 28 gauge with injecting into the i-port Advance™

Warnings and Precautions

- The i-port Advance™ is only sterile and non-pyrogenic if the insertion device is unopened or undamaged. DO NOT use if the package is already open or is damaged, or if the protective backing has

fallen off. Ensure the sterility by checking the sterile paper and tamper-proof seal.

- Carefully read through instructions prior to using the i-port Advance™ as failure to follow instructions may result in pain or injury.
- Inserter can be disposed of safely after insertion if the lid has been placed back on. For available sharps containers, please consult your local pharmacy.
- If i-port Advance™ is not securely placed in the inserter with the needle pointing straight ahead prior to insertion, pain or injury may occur.
- Be sure that the needle guard is removed before insertion.
- DO NOT use if introducer needle is bent or broken.
- DO NOT USE A NEEDLE LONGER THAN 5/16 in. or 8 mm (DO NOT attempt to use a 1/2 in. or 12 mm injection needle with the i-port Advance™. Doing so may damage the device causing unnecessary punctures of the skin, or tearing or puncturing of the soft cannula which may result in unpredictable medication delivery).
- DO NOT USE A NEEDLE SHORTER THAN 3/16 in. (5mm) when injecting into the i-port Advance™.
- The i-port Advance™ is a single use device. Reuse of the port may damage the cannula in the device and increases the risk of the adhesive tape becoming loose. Reuse of the port may result in infection and/or site irritation and inaccurate medication delivery.
- Inaccurate medication delivery, infection, and/or site irritation may result from improper application to the body and/or maintenance of the i-port Advance™ application site. If the application site becomes irritated or inflamed, remove and discard the device and apply a new device to a different location on the body.
- The device's soft cannula may crimp or bend when applying or while wearing. A crimped or bent cannula may limit or prevent medication delivery

through i-port Advance™. You should remove and discard the device immediately if you suspect a crimped or bent cannula. Indications that your i-port Advance™ has a crimped or bent cannula include:

- Difficulty pressing down on the plunger of the syringe when injecting
- Injected medication escaping from the device onto the top of the septum
- Injected medication leaking onto the skin or saturating the adhesive pad
- Moisture present within the clear body of the device
- DO NOT reinsert the introducer needle into the i-port Advance™. Reinsertion could cause puncturing or tearing of the soft cannula which would result in unpredictable medication delivery or injury.
- Replace the i-port Advance™ if the adhesive tape becomes loose or is displaced from original site. Since the cannula is soft, it will not cause pain if it slips out, and this may therefore take place without notice. The soft cannula must always be completely inserted to receive full amount of medication from the injection.
- Swab skin with alcohol or similar disinfectant where you intend to apply the i-port Advance™. (Let skin surface dry completely prior to applying the device to the disinfected area.)
- Ensure the septum is clean before every use (use pure water).
- DO NOT put perfumes or deodorants on i-port Advance™ as these may affect the integrity of the device.
- DO NOT inject into the i-port Advance™ through clothing.
- Never point loaded inserter towards any body part where insertion is not desired.

- Store i-port Advance™ at room temperature. DO NOT store or leave the device in direct sunlight.
- Avoid applying the i-port Advance™ to sites that contain an excess amount of scar tissue or lipoatrophy.
- DO NOT inject more than 75 times through a single device.
- DO NOT inject multiple drugs contraindicated for subcutaneous administration into the same site.
- The device has a residual volume up to 2.60uL (0.0026mL). It is not known if the residual drug remains in the cannula or is absorbed into the subcutaneous space over time. The prescribing healthcare provider should consider if the retention of, or delay in delivery of small doses of a drug (such as 1IU or 10uL of the U100 insulin) may adversely affect the therapeutic goals in some patients such as infants or small children.


Recommendations

- Test blood glucose per your healthcare provider's instructions to ensure that insulin delivery at the application site is appropriate.
- Test blood glucose 1 to 3 hours after introducing and using the i-port Advance™ to ensure insulin delivery performs as anticipated, and measure your blood glucose on a regular basis.
- Avoid applying and using a new i-port Advance™ prior to bedtime unless blood glucose can be checked 2 to 4 hours after injection.
- In the case of unexplainably high blood glucose, remove and discard the i-port Advance™ and apply a new device to a different location on the body since the soft cannula could be dislodged, crimped or partially clogged. Should suggested actions not provide a resolution to the problem or should your blood glucose remain high, contact your physician or healthcare provider.

APPLICATION

- 1 Wash hands with soap and water.
- 2 Swab skin with alcohol or similar disinfectant.
(Let skin surface dry completely prior to applying the device).
- 3 Pull red tab and remove seal.
- 4 Remove sterile paper.
- 5 Securely hold inserter and remove plastic lid.
- 6 Remove protective backing from adhesive pad.
- 7 Twist and remove needle guard from introducer needle.

Get to know the i-port *Advance*™

(A) Sterile Paper

(B) Inserter

(C) Injection Port

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(E) Plastic Lid

Center of Inserter (F)

Lined Ridges (G)

Round Indentations (H)

Introducer Needle (I)

Septum (J)


Adhesive Pad (K)

Soft Cannula (L)

- 8** With introducer needle facing down, put fingers on lined ridges then pull up on center of inserter until locked in upright position.
- 9** Place inserter against skin and put fingers on round indentations. Squeeze round indentations simultaneously to insert injection port.
- 10** Press down on center of inserter to secure injection port to skin.
- 11** Grab center of inserter only, and gently pull to remove inserter from injection port.
- 12** Secure adhesive pad to skin.
- 13** Replace plastic lid and dispose of according to local requirements.

USE

- 1** Clean septum before each use.
- 2** Puncture septum using injection needle.
- 3** Slowly inject medication.
- 4** Remove injection needle.

REMOVE

- 1** Grab edge of adhesive pad and remove from skin.

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