

Guardian™ app Official Compatibility Guide

Updated December 7, 2022



To use the Guardian™ 4 Smart CGM system, you will need a compatible device and operating system (OS), as listed below. It is important to check this compatibility table before installing a new OS or attempting to download the app on a new phone.

When a new OS is released, it may take time to complete testing. We recommend that you check this table before updating your OS.

Care partner SMS alerts are not limited to the below compatibility.

OUS	
Guardian™ 4 smart CGM system - Apple devices	
iOS	12.0, 12.1, 12.2, 12.3, 13.0, 13.1, 13.2, 13.3, 13.5.1, 13.6, 14.0, 14.2, 14.3, 14.4, 14.5, 14.6, 15.0, 15.1, 15.3, 15.4, 16
iPhone	6S Plus, SE (1 st Gen 2016), 7, 7 Plus, 8, 8 Plus, X, XS, XS Max, XR, 11, 11 Pro, 11 Pro Max, iPhone SE (2nd Gen 2020), iPhone 12, iPhone 12 mini, iPhone 12 Pro, iPhone 12 Pro Max, iPhone 13, iPhone 13 mini, iPhone 13 Pro, iPhone 13 Pro Max,
iPod Touch	7 th Generation
iPad	5,6, Air 2, Mini 4, Mini 5, Pro (9.7"), Pro (10.5"), Pro (12.9"), Pro 2 (12.9")

Guardian™ 4 smart CGM system - Android devices		
Brand	Model	Android OS
Samsung	S8 (SM-G950N, SM-G950F, SM-G950W, SM-G950U, SM-G9508, SM-G950U1, SM-G9500, SCV36, SC-02J)	8, 9
	S8+ (SM-G955U1, SM-G9550, SM-G955N, SM-G955F, SM-G955W, SM-G955U, SC-03J, SCV35)	
	S9 (SM-G960W, SM-G960U, SM-G960U1, SM-G960U, SCV38, SC-02K, SM-G9608, SM-G9600)	8, 9, 10

S9+ (SM-G965N, SM-G965F, SM-G965W, SM-G965U, SM-G965U1, SC-03K, SCV39)	
Note 8 (SC-01K, SCV37, SM-N950F, SM-N950N, SM-N950XN, SM-N950U, SM-N9500, SM-N9508, SM-N950W, SM-N950U1)	9
Galaxy A20 (SC-02M, SCV46, SCV46-J, SCV46-U, SM-A205F, SM-A205FN, SM-A205G, SM-205GN, SM-A205W, SM-A205YN, SM-A205U, SM-A205U1, SM-S205DL)	9, 10
Galaxy A11 (SM-A115A, SM-A115AP, SM-A115AZ, SM-A115F, SM-A115M, SM-A115U, SM-A115U1, SM-A115W)	10
Galaxy A21 (SC-42A, SCV49, SM-A215U, SM-A215U1, SM-A215W, SM-S215DL, SM-A217F, SM-A217M, SM-A217N)	10
Galaxy A50 (SM-A505F, SM-A505FM, SM-A505FN, SM-A505G, SM-A505GN, SM-A505GT, SM-A505N, SM-A505U, SM-A505U1, SM-A505W, SM-A505YN, SM-S506DL)	9, 10, 11
Galaxy A51 (SM-A515F, SM-A515U, SM-A515U1, SM-A515W, SM-S515DL) Galaxy A51 5G (SC-54A, SCG07, SM-A5160, SM-A516B, SM-A516N, SM-A516U, SM-A516U1, SM-A516V)	10, 11
Galaxy A71 (SM-A715F, SM-A715W) Galaxy A71 5G (SM-A7160, SM-A716B, SM-A716U, SM-A716U1, SM-A716V)	10, 11
S10 (SM-G973F, SM-G973N, SM-G9730, SM-G973U, SM-G973U1, SM-G973W, SC-03L, SCV41, SM-G9738, SM-G973C) S10 Lite (SM-G770F, SM-G770U1) S10e (SM-G970F, SM-G970N, SM-G9700, SM-G9708, SM-G970U, SM-G970U1, SM-G970W) S10+ (SM-G975F, SM-G975N, SM-G9750, SM-G975U, SM-G975U1, SM-G975W, SC-04L, SCV42, SM-G9758) S10+ Olympic Games Edition (SC-05L) S10 5G (SM-G977B, SM-G977N, SM-G977P, SM-G977T, SM-G977U)	9, 10, 11

	<p>S20 (SM-G980F)</p> <p>S20+ (SM-G985F)</p> <p>S20 5G (SC-51A, SCG01, SM-G9810, SM-G981N, SM-G981U, SM-G981U1, SM-G981V, SM-G981W, SM-G981B)</p> <p>S20+ 5G (SC-52A, SCG02, SM-G9860, SM-G986N, SM-G986U, SM-G986U1, SM-G986W, SM-G986B)</p> <p>S20 Ultra 5G (SCG03, SM-G9880, SM-G988N, SM-G988U, SM-G988U1, SM-G988W, SM-G988B)</p>	10, 11
	<p>S21 5G (SM-G9910, SM-G991U1, SM-G991U, SM-G991W, SM-G991B, SM-G991N)</p> <p>S21+ 5G (SM-G9960, SM-G996U1, SM-G996W, SM-G996B, SM-G996N)</p> <p>S21 Ultra 5G (SM-G9980, SM-G998U, SM-G998U1, SM-G998W, SM-G998B, SM-G998N)</p>	11,12
	<p>S22 5G (SM-S901B, SM-S901U, SM-S901U1, SM-S901W, SM-S901N, SM-S9010, SM-S901E)</p> <p>S22+ 5G (SM-S906B, SM-S906U, SM-S906U1, SM-S906W, SM-S906N, SM-S9060, SM-S906E)</p> <p>S22 Ultra 5G (SM-S908B, SM-S908U, SM-S908U1, SM-S908W, SM-S908N, SM-S9080, SM-S908E)</p>	12
	<p>Note 9 (SC-01L, SCV40, SM-N960F, SM-N960N, SM-N9600, SM-N960W, SM-N960U, SM-N960U1)</p>	10
	<p>Note 10 (SM-N970F, SM-N9700, SM-N970U, SM-N970U1, SM-N970W)</p> <p>Note 10 5G (SM-N971N)</p> <p>Note 10+ (SC-01M, SCV45, SM-N9750, SM-N975C, SM-N975U, SM-N975U1, SM-N975W, SM-N975F)</p> <p>Note 10+ 5G (SM-N976B, SSM-N976N, SM-N9760, SM-N976Q, SM-N976V, SM-N976U)</p>	10, 11
	<p>Note 20 (SM-N980F)</p>	10, 11

	Note 20 5G (SM-N9810, SM-N981N, SM-N981U, SM-N981U1, SM-N981W, SM-N981B) Note 20 Ultra (SM-N985) Note 20 Ultra 5G (SM-N9860, SM-N986N, SM-N986U, SM-N986U1, SM-N986W, SM-N986B)	
LG™	G8 ThinQ (LM-G820, LM-G820N) V50 ThinQ (LM-V450, LM-V500, LM-V500N) V50S ThinQ (LM-V510N) G8S ThinQ (LM-G810)	9, 10,
Xiaomi™	Mi 9, Mi MIX 2S	9, 10
Huawei™	P30 (ELE-L04, ELE-L09, ELE-L14, ELE-L29, ELE-L39, ELE-L49, ELE-TL00)	9, 10
Sony™	Xperia 1, Xperia 5, Xperia XZ2 (H8314, H8324)	9, 10, 11
Google	Pixel 3	9, 10, 11
Google	Pixel 3 XL, Pixel 3A, Pixel 3A XL, Pixel 4, Pixel 4A, Pixel 4A (5G), Pixel 4 XL, Pixel 5	10, 11, 12

*While the Guardian™ app is available for download on other models of Samsung, Xiaomi, Sony, LG and Huawei phones, it has not been fully tested and verified to be compatible.

